ottobock.

Quality for life

CE

C-Leg 3C98-3/3C88-3

DE	Gebrauchsanweisung (Benutzer)	3
EN	Instructions for use (user)	37
FR	Instructions d'utilisation (Utilisateur)	71
IT	Istruzioni per l'uso (Utilizzatore)	107
ES	Instrucciones de uso (Usuario)	143
PT	Manual de utilização (Usuário)	177
NL	Gebruiksaanwijzing (Gebruiker)	211
TR	Kullanma talimatı (Kullanıcı)	247

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Table of contents

EN

1	Foreword	39
2	Product Description	39
2.1	Design	39
2.2	Function	39
3	Application.....	40
3.1	Indications for use	40
3.2	Area of Application	40
3.3	Conditions of use	40
3.4	Indications	40
3.5	Qualification	40
4	Safety.....	41
4.1	Explanation of Warning Symbols	41
4.2	Structure of the safety instructions	41
4.3	General Safety Instructions.....	41
4.4	Information on the Power Supply/Battery Charging.....	43
4.5	Battery charger information.....	44
4.6	Information on Proximity to Certain Areas	44
4.7	Information on Use	45
4.8	Notes on the safety modes.....	46
4.9	Information on the use of a mobile device with the cockpit app.....	47
5	Scope of Delivery and Accessories.....	48
6	Charging the battery	48
6.1	Connecting the power supply and battery charger	48
6.2	Charging the prosthesis battery.....	49
6.3	Display of the current charge level	49
7	Cockpit app.....	50
7.1	System Requirements	50
7.2	Initial connection between cockpit app and prosthesis	51
7.2.1	Starting the cockpit app for the first time	51
7.3	Control elements for cockpit app	52
7.3.1	Cockpit app navigation menu.....	52
7.4	Managing prostheses.....	52
8	Use	53
8.1	Movement patterns in basic mode (mode 1)	53
8.1.1	Standing	53
8.1.1.1	Stance function	54
8.1.2	Walking.....	54
8.1.3	Sitting down	54
8.1.4	Sitting	55
8.1.5	Standing up	55

8.1.6	Walking up stairs.....	55
8.1.7	Walking down stairs	55
8.1.8	Walking down a ramp	56
8.1.9	Kneeling.....	56
8.2	MyModes	56
8.2.1	Switching MyModes with the cockpit app	56
8.2.2	Switching MyModes using motion patterns	57
8.2.3	Switching from a MyMode back to basic mode.....	58
8.3	Changing prosthesis settings.....	58
8.3.1	Changing the prosthesis setting using the cockpit app.....	59
8.3.2	Overview of adjustment parameters in basic mode.....	59
8.3.3	Overview of adjustment parameters in MyModes	60
8.4	Turning Bluetooth on the prosthesis on/off.....	60
8.4.1	Switching Bluetooth off/on using the cockpit app	60
8.5	Querying the prosthesis status	61
8.5.1	Query status through cockpit app	61
8.5.2	Status display in the cockpit app	61
8.5.3	Status information on the display of the remote control (optional accessory)	61
9	Additional operating states (modes)	61
9.1	Empty battery mode	61
9.2	Mode for charging the prosthesis	61
9.3	Safety mode.....	61
9.4	Overheating mode	62
10	Cleaning	62
11	Maintenance	62
12	Legal Information	62
12.1	Liability	62
12.2	Trademarks	62
12.3	CE Conformity	63
12.4	Local Legal Information	63
13	Appendices	64
13.1	Symbols Used.....	64
13.1.1	Symbols on the product	64
13.1.2	Symbols on the battery charger	65
13.2	Operating states/error signals.....	65
13.2.1	Signals for operating states	65
13.2.2	Warnings/error signals.....	66
13.2.3	Error messages while establishing a connection with the cockpit app	67
13.2.4	Status signals	68
13.3	Technical data.....	68

1 Foreword

INFORMATION

Last update: 2015-03-05

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Receive instruction in the correct and safe use of the product by qualified personnel.
- ▶ Please contact the qualified personnel if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The product "C-Leg" is called the product/prosthesis/knee joint/component in the following.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product Description

2.1 Design

The product consists of the following components:

1. Proximal pyramid adapter
2. LED (blue) as indicator for the Bluetooth connection
3. Battery and cover caps
4. Optional flexion stop (8°)
5. Hydraulic unit
6. Charging receptacle cover
7. Charging receptacle
8. Pylon stop
9. Distal tube clamp screws

2.2 Function

This product features microprocessor control of the stance and swing phase.

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the damping behaviour of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled stance and swing phase, the system can be individually adapted to your needs.

For this purpose, the product is adjusted by qualified personnel using adjustment software.

The product features MyModes for special motion types (e.g. cross-country skiing, ...). These are pre-configured by the prosthetist using the adjustment software and can be activated with special motion patterns, the cockpit app and the remote control (optional accessory) (see Page 56).

In case of a system malfunction, safety mode makes restricted operation possible. Predefined resistance parameters are configured in the product for this purpose (see Page 61).

Empty battery mode permits safe walking when the battery is drained. Resistance parameters that are predefined by the product are configured for this purpose (see Page 61).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

3 Application

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

3.2 Area of Application

Area of application according to the MOBIS mobility system:

The product is recommended for mobility grade 2 (restricted outdoor walker), mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands). Approved for a body weight of **136 kg max.**

3.3 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.).

Permissible ambient conditions are described in the technical data (see Page 68).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

3.4 Indications

- For users with knee disarticulation, transfemoral amputation and hip disarticulation
- For unilateral or bilateral amputation
- Dysmelia patients with residual limb characteristics corresponding to knee disarticulation, transfemoral amputation or hip disarticulation
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.5 Qualification

The product may be fitted only by qualified personnel authorised by Ottobock after completing the corresponding training.

4 Safety

4.1 Explanation of Warning Symbols

⚠ WARNING Warning regarding possible serious risks of accident or injury.

⚠ CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 Structure of the safety instructions

⚠ CAUTION

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General Safety Instructions

⚠ WARNING

Using the prosthesis while operating a vehicle

Accidents due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- All users are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.
- Observe national legal regulations for retrofitting your vehicle in accordance with the type of fitting.

INFORMATION

The ability of prosthesis wearers to operate a vehicle is determined on a case-by-case basis. Criteria include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, prosthesis design) and the amputee's individual abilities.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- Do not open the power supply unit, adapter plug or battery charger.
- Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ CAUTION

Failure to observe warnings/error signals

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- Observe the warnings/error signals (see Page 66) and corresponding change in damping settings.

⚠ CAUTION

Independent user manipulation of system components

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the prosthesis because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the prosthesis if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

⚠ CAUTION

Risk of pinching in the joint flexion area

Injuries due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the joint.

⚠ CAUTION

Penetration of dirt and humidity into the product

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that no solid particles or foreign objects can penetrate into the product.
- ▶ The knee joint is weather-proof but not resistant to corrosion. Therefore, the knee joint should not come into contact with salt water or chlorinated water. Do not use the knee joint under extreme conditions like diving or jumping into the water. The knee joint is not designed for permanent underwater use.
- ▶ After contact with water, remove the protector and hold the prosthesis with the sole of the foot facing up until the water has drained from the knee joint/tube adapter. Dry the knee joint and components with a lint-free cloth and allow the components to fully air dry.
- ▶ Should the knee joint come into contact with salt water or chlorinated water, it must be rinsed with fresh water immediately.
- ▶ The knee joint is not resistant against penetration from water jets or steam.
- ▶ If water gets into the electronics, the knee joint must be inspected by an authorised Ottobock Service Centre.

⚠ CAUTION

Mechanical stress during transport

- > Falling due to unexpected product behaviour as a result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- Only use the transport packaging for transportation.

⚠ CAUTION

Signs of wear on system components

Falling due to damage or malfunction of the product.

- In the interest of the user's own safety and in order to maintain the operational reliability and the validity of the warranty, the specified service intervals must be observed.

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock Dermaclean).

INFORMATION

When using exoprosthetic knee joints, servomotor, hydraulic, pneumatic or brake load dependent control functions can cause movement noise. This kind of noise is normal and unavoidable. It generally does not cause any problems. If movement noise increases noticeably during the lifecycle of the knee joint, the knee joint should be inspected by an authorised Ottobock Service Centre immediately.

4.4 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the product without taking it off

- > Falling due to walking and getting caught on a connected battery charger.
- > Falling due to unexpected product behaviour because of changed damping behaviour.
- For safety reasons, take the product off before charging the battery.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

⚠ CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- Check the power supply unit, charger and charger cable for damage before use.
- Replace any damaged power supply unit, charger or charger cable.

4.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see Page 68).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

4.6 Information on Proximity to Certain Areas

⚠ CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.50 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.35 m
 - DECT cordless phones incl. base station: 0.18 m
 - WiFi (routers, access points,...): 0.11 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.11 m

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
If this cannot be avoided, ensure at least that you are able to walk or stand steadily (e.g. by using a handrail or the support of another person).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected changes in the damping behaviour of the product.

CAUTION

Remaining in areas outside the allowable temperature range

Falling due to malfunction or the breakage of load-bearing product components.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see Page 68).

4.7 Information on Use

CAUTION

Walking up stairs

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking up stairs and place most of the area of the sole of your foot on the stair surface.
- ▶ Particular caution is required when carrying children up stairs.

CAUTION

Walking down stairs

Falling due to foot placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warnings/error signals (see Page 66).
- ▶ If pulsating vibration signals start (overheating mode activated), stop walking down stairs immediately and check whether stance phase control is active.
- ▶ Particular caution is required when carrying children down stairs.

CAUTION

Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- > Be sure to pay attention to pulsating vibration signals. They indicate the risk of overheating.
- > As soon as these pulsating vibrations begin, you have to reduce the activity level so the hydraulic unit can cool down.
- > Full activity may be resumed after the pulsating vibration signals stop.
- > If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. The product should be inspected by an authorised Ottobock Service Centre.

⚠ CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, paragliding, etc.).
- Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by a prosthetist. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- Ensure that you stand securely during all switching processes.
- Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- Switch back to basic mode once the activities in MyMode have been completed.
- Take the weight off the product and correct the switching, if required.

⚠ CAUTION

Improper use of the stance function

Falling due to unexpected product behaviour because of changed damping behaviour.

- Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the prosthesis.
- Make sure you have been instructed in the proper use of the stance function by authorised Ottobock personnel. Information on the stance function see Page 54.

4.8 Notes on the safety modes

⚠ CAUTION

Using the product in safety mode

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- Observe the warnings/error signals (see Page 66).

⚠ CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- Do not continue using the defective product.
- Contact your prosthetist immediately.

⚠ CAUTION

Safety mode cannot be deactivated

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If you cannot deactivate safety mode by recharging the battery, a permanent malfunction has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.

⚠ CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warning/error signals (see Page 66).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. The prosthetist is your contact person.

4.9 Information on the use of a mobile device with the cockpit app

⚠ CAUTION

Improper use of the device

Falling due to altered damping behaviour as a result of unexpected switching into MyMode.

- ▶ Make sure you have been instructed in the proper use of the device with the cockpit app.

⚠ CAUTION

Independently applied changes or modifications made to the device

Falling due to altered damping behaviour as a result of unexpected switching into MyMode.

- ▶ Do not make any independent changes to the hardware of the device.
- ▶ Do not make any independent changes to the software/firmware of the device which are not included in the update function of the software/firmware.

⚠ CAUTION

Improper mode switching with the device

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter (beeper) and the display on the device.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.

NOTICE

Destruction of the device due to falling or penetration of water

Malfunction of the device.

- ▶ If necessary, let the device dry at room temperature (at least 1 day).
- ▶ If you can no longer switch from MyMode back to basic mode, the only way to switch the component back to basic mode is by using a motion pattern (see Page 58) or plugging/unplugging the battery charger.

NOTICE

Failure to observe the system requirements for the installation of the cockpit app

Malfunction of the device.

- ▶ Install the cockpit app only on the operating systems listed in the section "System requirements" (see Page 50). The tested devices are listed in this section as well.

5 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. 3C88-3 C-Leg knee joint (with threaded connector) or
- 1 pc. 3C98-3 C-Leg (with pyramid connector)
- 1 pc. 757L16* power supply
- 1 pc. 4E50 Battery Charger for C-Leg
- 1 pc. cosmetic case for battery charger and power supply
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. 646D790, 646D790=1 instructions for use (user)
- 1 pc. "Cockpit 4X441-Andr=V*" Android app

Accessories

The following components are not included in the scope of delivery and may be ordered separately:

- 3S26 cosmetic foam cover
- 4X860=* C-Leg Protective Cover (w/o shield)
- 4P863* Shield Insert
- 4X156 Charger Extension Cable – Ankle
- 4X157 Charger Extension Cable – Knee
- 4X350 Remote Control

6 Charging the battery

The following points must be observed when charging the battery:

- With uninterrupted walking, the capacity of the fully charged battery is sufficient for at least 16 hours. It lasts about 2 days with average use.
- We recommend charging the product overnight when used on a daily basis.
- When used daily, the complete charging unit (power supply – battery charger) may remain plugged into the wall socket.
- The battery should be charged for at least 3 hours prior to initial use.
- Use the 757L16* power supply and 4E50* battery charger to charge the battery.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place.

- 2) Connect the round, **three-pin** plug of the power supply to the **12V** receptacle on the battery charger so that the plug locks into place.
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 3) Connect the round, **four-pin** plug of the charging cable to the **OUT** receptacle on the battery charger so that the plug locks into place.
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 4) Plug the power supply unit into the wall socket.
 - The green LED on the back of the power supply and the green LED on the battery charger light up.
 - If the green LED on the power supply and the green LED on the battery charger do not light up, there is an error (see Page 66).

6.2 Charging the prosthesis battery

- 1) Open the charging receptacle cover.
- 2) Connect the charging plug to the charging receptacle of the product.
INFORMATION: Take note of the plug direction!
 - The correct connection of the battery charger to the product is indicated by feedback (see Page 68).
- 3) The charging process starts.
 - Once the product battery is fully charged, the yellow LED on the battery charger turns off.
- 4) Disconnect the battery charger after the charging process is complete.
 - A self-test is performed. The joint is operational only after corresponding feedback (see Page 68).
- 5) Close the charging receptacle cover.

6.3 Display of the current charge level

INFORMATION

The charge level cannot be displayed during the charging process.

- 1) Rotate the prosthesis 180° (the sole of the foot must face up).
- 2) Hold still for 2 seconds and wait for beeps.

Beep signal	Vibration signal	Battery charge level
5x short		more than 80%
4x short		65% to 80%
3x short		50% to 65%

Beep signal	Vibration signal	Battery charge level
2x short		35% to 50%
1x short	3x long	20% to 35%
1x short	5x long	less than 20%

Display of the current charge level using the cockpit app:


Once the cockpit app has been started, the current charge level is displayed in the bottom line of the screen:

1. 38% – battery charge level of currently connected prosthesis

Display of current charge level using remote control (optional accessory):

The current charge level is displayed in the status bar of the remote control when it is connected and switched on:

1. – Charge level of battery for currently connected prosthesis
2. – Prosthesis is being charged
2. 68% – Charge level of battery for currently connected prosthesis, as a percentage

7 Cockpit app

The cockpit app enables switching from basic mode into the pre-configured MyModes. In addition, information about the product (step counter, battery charge level, etc.) can be called up.

The behaviour of the product can be changed to a certain extent on a day-to-day basis using the app (e.g. while becoming accustomed to the product). At the next appointment, the prosthetist can trace the changes using the adjustment software.

INFORMATION

The cockpit app can be downloaded free of charge from the respective online store. For more information, please visit the website: <http://www.ottobock-group.com/cockpitapp>. You can also read the QR code on the supplied Bluetooth PIN card.

7.1 System Requirements

The functioning of the cockpit app is assured on devices with the Android operating system version 4.0.3 or higher.

The functionality was tested on the following devices:

- Samsung Galaxy S5, Galaxy S4, Galaxy S4 mini, Galaxy SIII, Galaxy SIII mini, Galaxy Note II, Galaxy Fame
- Sony Xperia Z, Xperia Z3, Xperia J, Xperia SP
- HTC One, One mini
- LG Optimus L9, Optimus F5, Optimus F6, Optimus G, Optimus G2
- Huawei Ascend P6, Ascend G500
- Motorola Droid Razr Maxx, Moto X, Nexus 6

7.2 Initial connection between cockpit app and prosthesis

The following points should be observed before the initial connection:

- Bluetooth of the component must be switched on.
If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down (sole of the foot must point up) or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 60).
- Bluetooth of the device must be switched on.
- The device must not be in "flight mode" (offline mode), otherwise all wireless connections are turned off.
- The device must be connected to the internet.
- The serial number and PIN code of the component being connected must be known. They are found on the enclosed card. The serial number begins with the letters "SN".


7.2.1 Starting the cockpit app for the first time

- 1) Tap the  symbol.
→ The end user license agreement (EULA) is displayed.
 - 2) Accept the end user license agreement (EULA) by tapping the **Accept** button. If the end user license agreement (EULA) is not accepted, the cockpit app cannot be used.
→ The welcome screen appears.
 - 3) Tap the **Add component** button.
→ The "**Preparation**" screen appears. The serial number of the component should be entered here.
 - 4) Follow the subsequent instructions on the screen.
 - 5) After the PIN code is entered, a connection to the prosthesis is established.
→ While the connection is being established, 3 beep signals sound and the  symbol appears.
The  symbol is displayed when the connection has been established.
- Once the connection has been established successfully, the data are read from the prosthesis. This process may take up to a minute.
The main menu will then appear with the name of the connected prosthesis.

INFORMATION


After the initial connection to the prosthesis has been established successfully, the app will connect automatically each time it is started. No further steps are required.

7.3 Control elements for cockpit app

1. ☰ Call up the navigation menu (see Page 52)
2. Product
The component name can only be changed with the adjustment software.
3. If connections to more than one component have been saved, you can switch between the saved components by tapping the entry **change**.
4. MyModes configured with the adjustment software.
Switching the mode by tapping the corresponding symbol and confirming by tapping "**OK**".
5. Currently selected mode
6. Charge level of the component.
 - Component battery fully charged
 - Component battery emptyThe current charge level is also displayed in %.
7. ⚡ Display and name of the currently selected mode (e.g. **1. Basic Mode**)
8. (i:i) Connection to component has been established
(○) Connection to component has been interrupted. The app is attempting to re-establish the connection automatically.

7.3.1 Cockpit app navigation menu

Tap the ☰ symbol in the menus to display the navigation menu. Additional settings for the connected component can be configured in this menu.

Product

Name of the connected component

MyModes

Return to the main menu to switch MyModes

Functions

Call up additional functions of the component (e.g. turn off Bluetooth) (see Page 60)

Settings

Change settings of the currently selected mode (see Page 58)

Status

Query status of the connected component (see Page 61)

Manage components

Add or delete components (see Page 52)

Imprint/Info

Display information/legal notices for the cockpit app

7.4 Managing prostheses

Connections with up to four different components can be stored in the app. However, a component can only be connected to one device or remote control at a time.

Adding component

- 1) Tap the ☰ symbol in the main menu.
→ The navigation menu opens.

- 2) In the navigation menu, tap the entry "**Manage components**".
 - 3) In the following screen, tap the "**Add component**" button.
 - The "Preparation" screen appears. The serial number of the component should be entered here. The serial number begins with the letters "SN".
 - 4) Follow the subsequent instructions on the screen.
 - 5) After the PIN code is entered, a connection to the prosthesis is established.
 - While the connection is being established, 3 beep signals sound and the  symbol appears.
 - The  symbol is displayed when the connection has been established.
 - Once the connection has been established successfully, the data are read from the prosthesis. This process may take up to a minute.
- The main menu will then appear with the name of the connected prosthesis.


Deleting a component

- 1) Tap the  symbol in the main menu.
 - The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Manage components**".
- 3) Then tap the option "**Manage components**".
- 4) Tap the  symbol under the component you want to delete.
 - The component is deleted.

Connecting component with multiple devices

A component can be connected to multiple devices (Managing prostheses).

If there is an existing connection between the component and a different device, the following information appears while the connection is being established:


- ▶ Tap the **OK** button.
- The connection to the last connected device is broken off and established with the current device.

8 Use

8.1 Movement patterns in basic mode (mode 1)

8.1.1 Standing

Knee control through high hydraulic resistance and static alignment.

A stance function can be enabled using the adjustment software. Please see the following section for further information on the stance function.

8.1.1.1 Stance function

INFORMATION

If this function has been enabled in the adjustment software, the stance function can be activated/deactivated using the cockpit app (see Page 59) or the remote control (optional accessory).

The stance function (standing mode) is a functional supplement to the basic mode (mode 1). This function makes it easier for the user to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction at a flexion angle between 5° and 65°.

This function must be enabled in the adjustment software. Once the function is enabled, you can also choose between an automatic and manual lock.

Auto joint lock

The auto stance function recognises any situation that puts strain on the prosthesis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The knee joint is always locked in the flexion direction when the prosthetic leg is not fully extended and is at rest. Upon forward or backward rollover or extension, the level of resistance is immediately reduced to stance phase resistance again.

The knee joint is not locked when the above conditions are met and a sitting position is assumed (for example while driving).

The stance function can also be set manually to any flexion angle between 5° und 65° (see the next section).


Manual joint lock

- 1) Flex the joint between 5° and 65° and keep it still for one second.
- 2) Slowly extend the joint up to the desired standing angle.
- 3) In this position, keep the joint still for one second until it vibrates.
→ The blocked joint can now be fully loaded in the flexion direction.

Releasing joint lock

- Quickly extend the joint or tilt the thigh slightly forward or more than 50° back.

8.1.2 Walking

Initial attempts at walking with the prosthesis always require the instruction of trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires that the prosthesis roll over to the front out of the stride position.

8.1.3 Sitting down

The resistance of the product while sitting down ensures that the knees bend evenly, thereby supporting the contralateral side.

- 1) Position both feet side by side at the same level.
- 2) While sitting down, the user should distribute his/her weight evenly between both legs and use armrests, if available.
- 3) Move the buttocks in the direction of the backrest and lean the upper body forward.

8.1.4 Sitting

INFORMATION

While sitting, the knee joint switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is enabled in the adjustment software or not.

If the patient is in a sitting position for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the knee joint switches the resistance to a minimum in the extension direction.

If the sitting function is activated in the adjustment software and switched on via the cockpit app or remote control (optional accessory) (see Page 59), resistance is also reduced in the flexion direction.

8.1.5 Standing up

Flexion resistance is increased steadily while standing up.

- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Put the hands on armrests, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on feet.


8.1.6 Walking up stairs

Walking up stairs step-over-step is not possible.

- 1) Hold the handrail with one hand.
- 2) Place the foot of the contralateral leg on the first step.
Bring up the leg with the prosthesis.


8.1.7 Walking down stairs

This function must be practised and executed consciously. Only when the sole is properly positioned can the system react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the product on the step so that the foot projects halfway over the edge of the step.
→ This is the only way to assure a secure rollover.
- 3) Roll the foot over the edge of the step.
→ This flexes the product slowly and evenly at the knee joint.
- 4) Place the foot of the contralateral leg on the next step.
- 5) Place the foot of the prosthesis leg on the next step after that.

8.1.8 Walking down a ramp

Under increased flexion resistance, permit controlled flexion of the knee joint which lowers the body's centre of gravity.
The swing phase is not triggered even though the knee joint is flexed.



8.1.9 Kneeling

Under increased flexion resistance, permit controlled flexion of the knee joint to gradually reach the kneeling position. A hard impact of the knee joint on the ground should be avoided so the electronics are not damaged.
Using the 4X860= C-Leg protector is recommended for kneeling frequently.

8.2 MyModes

With the help of adjustment software, the prosthetist can activate and configure MyModes in addition to the basic mode. These can be selected by using the cockpit app, the remote control (optional accessory) or motion patterns. Switching by using motion patterns has to be activated in the adjustment software by the prosthetist.


These modes are intended for specific motion patterns or postures (e.g. inline skating, ...). Settings can be adjusted using the cockpit app or remote control (see Page 60).

8.2.1 Switching MyModes with the cockpit app

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the cockpit app.
If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 60).

Once a connection to a prosthesis has been established, the cockpit app can be used to switch between the MyModes.

- 1) Tap the symbol of the MyMode (1) you want in the main menu of the app.
→ A security question for changing the MyMode appears.
- 2) If you want to change the mode, tap the "OK" button.
→ A beep signal sounds to confirm the switch.
- 3) After switching, a symbol (2) is displayed to identify the active mode.
→ The current mode is also indicated by the symbol and name on the lower edge of the screen (3).

8.2.2 Switching MyModes using motion patterns

Information on switching

- Switching and the number of motion patterns must be activated by the prosthetist in the adjustment software.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Switching using motion patterns has to be enabled in the adjustment software.
- Position the prosthetic leg back slightly (lunge position) and bounce on the forefoot while maintaining constant contact with the floor.
Then keep the prosthetic leg still in this position (lunge position) for about 1 second without lifting the leg. Taking the weight off is no longer necessary.

Switching process

- 1) Position the prosthetic leg back slightly (lunge position).
- 2) While maintaining constant contact with the floor, bounce on the forefoot a number of times in one second depending on the desired MyMode (MyMode 1 = 3 times, MyMode 2 = 4 times).
- 3) **Keep the prosthetic leg still in this position (lunge position) for about 1 second without lifting the leg.** Taking the weight off is no longer necessary.
→ A beep signal will sound to confirm that the motion pattern has been recognised.
INFORMATION: If this beep signal does not sound, the requirements were not observed when bouncing the foot.
- 4) Move the prosthetic leg in next to the contralateral leg, set it down and keep it still for approx. 1 second.
→ A confirmation signal will sound to indicate that the prosthesis has successfully switched to the corresponding MyMode (2 times = MyMode 1, 3 times = MyMode 2).
INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly repositioned and kept still. Repeat the process to correctly switch to the required mode.

8.2.3 Switching from a MyMode back to basic mode

Information on switching

- Regardless of the configuration of additional MyModes in the adjustment software, it is always possible to switch back to basic mode (mode 1) with a motion pattern.
- It is always possible to switch back to basic mode (mode 1) by connecting/disconnecting the battery charger.
- Before the first step, always check whether the selected mode corresponds to the required motion type.

Requirements for successful switching using motion patterns

The following points must be observed to carry out switching successfully:

- Position the prosthetic leg back slightly (lunge position) and bounce on the forefoot while maintaining constant contact with the floor.

Then keep the prosthetic leg still in this position (lunge position) for about 1 second without lifting the leg. Taking the weight off is no longer necessary.

Switching process

- 1) Position the prosthetic leg back slightly (lunge position).
- 2) While maintaining constant contact with the floor, bounce on the forefoot at least 3 times but not more than 5 times.
- 3) **Keep the prosthetic leg still in this position (lunge position) for about 1 second without lifting the leg.** Taking the weight off is no longer necessary.
→ A beep signal will sound to confirm that the motion pattern has been recognised.

INFORMATION: If this beep signal does not sound, the requirements were not observed when bouncing the foot.

- 4) Move the prosthetic leg in next to the contralateral leg, set it down and keep it still for approx. 1 second.

→ A confirmation signal will sound to indicate that the prosthesis has successfully switched over to basic mode.

INFORMATION: If this confirmation signal does not sound, the leg with the prosthesis was not correctly repositioned and kept still. Repeat the process to correctly switch to the required mode.

8.3 Changing prosthesis settings

Once an active connection to a prosthesis has been established, the settings of the mode currently active can be changed using the cockpit app or the remote control (optional accessory).

INFORMATION

Bluetooth on the prosthesis must be switched on to change the prosthesis settings.


If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. The connection must be established during this period.

Information for changing the prosthesis settings

- Before changing settings, always check the main menu of the cockpit app or the display of the remote control (optional accessory) to make sure the correct prosthesis has been selected. Otherwise, parameters could be changed for the wrong prosthesis.
- It is not possible to change prosthesis settings nor to switch to a different mode while the prosthesis battery is being charged. Only the status of the prosthesis can be called up. Instead of the  symbol, the  symbol appears in the bottom row of the screen in the cockpit app.

- The prosthetist's setting is located in the middle of the scale. After making adjustments, this setting can be restored by tapping the "**Standard**" button (cockpit app) or by moving the slider control into the middle (remote control).
- Prostheses settings should be optimised using the adjustment software. The cockpit app or remote control (optional accessory) is not intended for use by the prosthetist in order to set up the prosthesis. The behaviour of the prosthesis can be changed to a certain extent on a day-to-day basis using the remote control or app (e.g. while becoming accustomed to the prosthesis). The prosthetist can use the adjustment software to trace these changes at the next appointment.
- If the settings of a MyMode are to be modified, one must first switch to this MyMode.

8.3.1 Changing the prosthesis setting using the cockpit app

- Once the component is connected and in the desired mode, tap the **E** symbol in the main menu.
→ The navigation menu opens.
 - Tap the menu entry "**Settings**".
→ A list appears with the parameters for the currently selected mode.
 - Change the setting of the desired parameter by tapping the "<", ">" symbols.
- INFORMATION:** The prosthetist's setting is marked and, after the setting has been changed, can be restored by tapping the "**Standard**" button.

8.3.2 Overview of adjustment parameters in basic mode

The parameters in basic mode describe the dynamic behaviour of the prosthesis in a normal gait cycle. These parameters act as basic settings for automatically adjusting the damping behaviour to the current motion situation (e.g. ramps, slow walking speed, etc.).

The stance function and/or the sitting function can also be activated/deactivated. Further information on the stance function (see Page 54). Further information on the sitting function (see Page 55).

The following parameters can be modified:

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Resistance	120 to 180	+/- 10	Resistance against flexion motion e.g. when walking down stairs or when sitting down

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Stance function		0 - deactivated 1 - activated	Activation/deactivation of the stance function. This function needs to be enabled in the adjustment software.
Sitting function		0 - deactivated 1 - activated	Activation/deactivation of the sitting function. This function needs to be enabled in the adjustment software.
Volume	0 to 4	0 to 4	Volume of beep signal for confirmation tones.

8.3.3 Overview of adjustment parameters in MyModes

The parameters in the MyModes describe the static behaviour of the prosthesis for a specific motion pattern such as cross-country skiing. Damping behaviour is not automatically controlled and adjusted in MyModes.

The following parameters can be modified in MyModes:

Parameter	Adjustment software range	App/remote control adjustment range	Meaning
Gain	0 to 100	+/- 10	Value for speed at which flexion resistance increases as the knee angle increases
Volume	0 – 4	0 – 4	Volume of beep signal for confirmation tones.

8.4 Turning Bluetooth on the prosthesis on/off

INFORMATION

Bluetooth on the prosthesis must be turned on in order to use the cockpit app. If Bluetooth is switched off, it can be turned on by turning the prosthesis upside-down or by connecting/disconnecting the battery charger. Bluetooth is then turned on for approx. 2 minutes. During this time, the app must be started and used to establish a connection. If required, Bluetooth on the prosthesis can be switched on permanently afterwards (see Page 60).

8.4.1 Switching Bluetooth off/on using the cockpit app

Switching off Bluetooth

- 1) When the component is connected, tap the  symbol in the main menu.
→ The navigation menu opens.
- 2) In the navigation menu, tap the entry "**Functions**".
- 3) Tap the entry "**Deactivate Bluetooth**".
- 4) Follow the on-screen instructions.

Switching on Bluetooth

- 1) Rotate the component or connect/disconnect the battery charger.
→ Bluetooth is switched on for approx. 2 minutes. During this time, the app must be started in order to establish a connection to the component.
- 2) Follow the on-screen instructions.

→ If Bluetooth is switched on, the  symbol appears on the screen.

8.5 Querying the prosthesis status

8.5.1 Query status through cockpit app

- 1) When the component is connected, tap the  symbol in the main menu.
- 2) In the navigation menu, tap the entry "Status".

8.5.2 Status display in the cockpit app

Menu option	Description	Possible actions
Trip: 1747	Daily step counter	Reset the counter by tapping the "Reset" button.
Step: 1747	Total step counter	For informational purposes only
Batt.: 68	Current prosthesis charge level, as a percentage	For informational purposes only

8.5.3 Status information on the display of the remote control (optional accessory)

Menu option	Description	Possible actions
Trip:1747	Daily step counter	Reset the counter by using the  key to confirm the menu option
Step:1747	Total step counter	For informational purposes only
Batt.:68	Current prosthesis charge level, as a percentage	For informational purposes only
Bluetooth: On	Switch Bluetooth function of the prosthesis on or off	The Bluetooth function of the prosthesis can be turned on or off (see Page 60) by confirming the menu item with the  key.

9 Additional operating states (modes)

9.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level is 0% (see Page 66). During this time, damping settings are set to their safety mode values. The prosthesis is then switched off. You can switch back to basic mode (mode 1) from empty battery mode by charging the product.

9.2 Mode for charging the prosthesis

The product is non-functional during charging.

The product is set to safety mode damping. This may be low or high depending on the setting in the adjustment software.

9.3 Safety mode

The product automatically switches to safety mode if a critical system fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

Default damping values are activated in safety mode. This makes limited walking possible for the user even though the system is not active.

The switch to safety mode is indicated by beeps and vibration signals immediately prior to switching (see Page 66).

Safety mode can be disabled by connecting then disconnecting the battery charger. If the product switches into safety mode again, this means a permanent error exists. The product must be inspected by an authorised Ottobock Service Centre.

9.4 Overheating mode

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), damping is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the damping settings that existed before the overheating mode.

Overheating mode is indicated by a long vibration every 5 seconds.

When overheating mode is active, it is not possible to switch to a MyMode, display the current battery charge level or change the prosthesis settings. In addition, the resistance in the flexion direction is not reduced when the sitting function is active.

Overheating mode is not activated in the MyModes.

10 Cleaning

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

11 Maintenance

INFORMATION

This component was tested for three million load cycles in accordance with ISO 10328.

Depending on the activity level, this corresponds to a service life of three to five years.

The duration of use can be individually extended depending on the intensity of use by performing regular service inspections.

Regular service inspections are recommended in the interest of the user's own safety, and in order to maintain operating reliability and protect the warranty. These service inspections include an inspection of the sensors and replacement of worn parts.

To have a service inspection carried out, please send the product as well as the battery charger and power supply unit to an authorised Ottobock Service Centre.

12 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

12.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

12.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

This product meets the requirements of the European Directive 1999/5/EC for radio equipment and telecommunications terminal equipment. The conformity assessment was drawn up by the manufacturer in accordance with Annex IV of the directive.

12.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA
Phone + 1-801-956-2400
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

13 Appendices

13.1 Symbols Used

13.1.1 Symbols on the product

Legal manufacturer

Compliance with the requirements according to "FCC Part 15" (USA)

Compliance with the requirements according to the "Radiocommunications Act" (AUS)

Non-ionising radiation

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN YYYY WW NNN Serial number

IP67

Dustproof, protection against temporary submersion

Caution, hot surface

13.1.2 Symbols on the battery charger

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

LOT PPPP YYYY WW Lot number

13.2 Operating states/error signals

The prosthesis indicates operating states and error messages through beeps and vibration signals.

13.2.1 Signals for operating states

Battery charger connected/disconnected

Beep signal	Vibration signal	Event
1x short		Battery charger is connected or battery charger already disconnected prior to start of charging mode
	3x short	Charging mode started (3 sec. after connection of battery charger)
1x short	1x before beep signal	Battery charger disconnected after start of charging mode

Mode switching

Beep signal	Vibration signal	Additional action performed	Event
1x short	1x short	Mode switching using remote control	Mode switching process using remote control.
1x short	1x short	Bouncing on the forefoot	Bouncing pattern recognised.
1x short	1x short	Prosthetic leg moved next to contralateral leg, set down and kept still for 1 second	Switching to basic mode (mode 1) carried out.
2x short	2x short	Prosthetic leg moved next to contralateral leg, set down and kept still for 1 second	Switching to MyMode 1 (mode 2) carried out.

Beep signal	Vibration signal	Additional action performed	Event
3x short	3x short	Prosthetic leg moved next to contralateral leg, set down and kept still for 1 second	Switching to MyMode 2 (mode 3) carried out.

13.2.2 Warnings/error signals

Error during use

Beep signal	Vibration signal	Event	Required action
	1 x long at interval of approx. 5 seconds	Overheated hydraulic unit	Reduce activity.
	3x long	Battery charge level under 25%	Charge battery soon.
	5x long	Battery charge level under 15%	Charge battery promptly, since the product is switched off after the next warning signal.
10x long	10x long	Battery charge level 0% After the beep and vibration signal, the product switches to empty battery mode and then switches off.	Charge the battery.
30x long	1x long, 1x short repeated every 2.6 seconds	Severe error e.g. failure of valve drives Possibly no switching into safety mode.	Restricted walking possible. Attempt to reset this error by connecting/disconnecting the battery charger. The battery charger must be connected for at least 5 seconds before it is disconnected. If the error persists, it is no longer permissible to use the product. The product must be inspected immediately by a prosthetist.

Error while charging the product

LED on power supply	LED on battery charger	Error	Resolution
		Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.

LED on power supply	LED on battery charger	Error	Resolution
○	○ ○ ⓘ	Non-functional socket	Check socket with another electric appliance.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	○ ○ ⓘ	No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
●	○ ● ⓘ	Battery is fully charged (or connection with product is interrupted).	<p>Take note of the beep signals to differentiate between these states.</p> <p>When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a single beep and vibration signal.</p> <p>The battery is fully charged if this beep signal is heard.</p> <p>If no beep signal is heard, the connection to the product is interrupted.</p>
			<p>If the connection to the product is interrupted, an authorised Ottobock Service Centre must inspect the product, battery charger and power supply.</p>

13.2.3 Error messages while establishing a connection with the cockpit app

Error message	Cause	Correction
Device not found	The connection could not be established as no component with the entered serial number could be found.	Compare the serial number entered with that of the component, then try to establish the connection again.
Component was connected to another device. Establish connection?	The component was connected to another device/remote control	To disconnect the original connection, tap the "OK" button. If the original connection is not to be disconnected, tap the "Cancel" button.

Error message	Cause	Correction
	A current connection to the prosthesis was interrupted	<p>Check the following points:</p> <ul style="list-style-type: none"> Distance from the prosthesis to the device Battery charge level of the prosthesis Bluetooth on the prosthesis switched on? (see Page 60) If multiple prostheses were stored in the remote control, was the correct prosthesis selected?

13.2.4 Status signals

Battery charger is connected

LED on power supply	LED on battery charger	Event
		Power supply and battery charger operational

Battery charger disconnected

Beep signal	Vibration signal	Event
1x short	1x short	Self-test completed successfully. Product is operational.
3x short		<p>Maintenance note Conduct the self-test again by connecting/disconnecting the battery charger. If the beep signal sounds again, visit your prosthetist soon. If necessary, the prosthetist will forward the product to an authorised Ottobock Service Centre.</p> <p>The product can be used without restrictions. However, vibration signals may not be generated.</p>

Battery charge level

Battery charger	
	Battery is being charged, battery charge level is less than 50%
	Battery is being charged, battery charge level is over 50%
	<p>Battery is fully charged (or connection with product is interrupted). Take note of the beep signals to differentiate between these states. When the battery charger is connected or disconnected, a self-test is conducted and confirmed by a single beep and vibration signal.</p> <p>The battery is fully charged if this beep signal is heard. If no beep signal is heard, the connection to the product is interrupted.</p>

13.3 Technical data

Environmental conditions	
Storage and transport in original packaging (≤ 3 months)	-20°C/-4°F to +40°C/+104°F

Environmental conditions	
Storage and transport without packaging (<48 hours)	-25°C/-13°F to +70°C/+158°F Max. 93% relative humidity, non-condensing
Long-term storage (>3 months)	-20°C/-4°F to +20°C/+68°F Max. 93% relative humidity, non-condensing
Operation	-10°C/+14°F to +60°C/+140°F Max. 93% relative humidity, non-condensing
Charging the battery	+10°C/+50°F to +45°C/+113°F

Product	
Article number	3C98-3*/3C88-3*
Mobility grade according to MOBIS	2 to 4
Maximum body weight	136 kg/300 lb
Protection rating	IP67
Range of Bluetooth connection to PC	max. 10 m/32 ft
Range of Bluetooth connection to remote control/mobile device	max. 10 m/32 ft
Maximum possible flexion angle	130°
Maximum possible flexion angle with pre-installed flexion stops	122°
Weight of the prosthesis without Protector	approx. 1235 g/43.56 oz

Prosthesis battery	
Battery type	Li-Ion
Charging cycles (charging and discharging cycles) after which at least 80% of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is non-functional
Operating time of the prosthesis with new, fully charged battery at room temperature	At least 16 hours of uninterrupted walking Approx. 2 days with average use

Battery charger	
Reference number	4E50*
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F max. 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Input voltage	12 V ==

Cockpit app	
Article number	Cockpit 4X441-Andr=V*

Cockpit app	
Supported operating system	Android 4.0.3 or higher
Website for download	http://www.ottobock-group.com/cockpitapp

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

— Consult the dealer or an experienced radio/ TV technician for help.
Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.
5470 Harvester Road
L7L 5N5 Burlington, Ontario
Canada
Phone + 1-800-665-3327

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Otto Bock Healthcare Products GmbH